

SN200 Non Surgical Isolation Gowns

AAMI Level 1, Level 2, Level 3

SN200 Non Surgical Isolation Gowns

AAMI LEVEL 1, LEVEL 2, LEVEL 3

Gowns are used to protect the wearer from the spread of infection or illness if the wearer comes in contact with potentially infectious liquid and solid material. Gowns are one part of an overall infection-control strategy.

Super Safe provides, among several apparel solutions, high quality AAMI Level 1, AAMI Level 2 and AAMI Level 3 full back isolation gowns, made of highly efficient filtering protective materials. The Super Safe SN200 Infection control gowns are comfortable, lightweight and breathable. Tear-resistant, flame retardant, and anti-fluid.

The SN200 gowns will protect staff members from blood, urine, saline, or other chemicals and body fluids during surgical procedures.

AAMI LEVEL 2

- Water resistance: Impact penetration AATCC 42
- · Water resistance: Hydrostatic pressure

AATCC 42

- Measures the resistance of fabrics to the liquid penetration of water by impact.
- Spray impact ≤ 1.0g
- Hydrostatic pressure ≥ 20cm
- · Low fluid levels

Product Details	
Product name	Super Safe SN200
Disinfecting Type	AAMI Level 1, AAMI Level 2, AAMI Level 3
Stock	Yes
Quality Certification	ISO13485/EN13795/SGS test report
Sterile:	Non – sterile
Application:	Safety protection
Place of Origin:	Israel
Туре	Standard, reinforced
Size	M-S/M/L/XL/XXL
Shelf Life	3 years
Material	SMS
Color	Blue

Materials descriptions

The Super Safe SN200 Infection control isolation gowns are produced out of strong and breathable SMS material (Spunbond/Meltblown/Spunbond), a multi-layer fabric composed of inner layers of meltblown polypropylene between outer layers of spunbond polypropylene that is ideal for extended wear. Light/medium SMS is appropriate for low amounts of fluid, while heavy weight SMS is more appropriate for moderate contact with fluids.

Features

- Full back disposable surgical gown
- · Single-use
- Secure protection
- Anti-fluid, anti-alcohol, anti-bloodm, anti-static
- Durable
- Comfortable, lightweight and breathable
- · Tear-resistant
- Flame retardant

Product applications Style guide



Product codes according to size				
Size	S (6.5 – 7)	M (7.5 – 8)	L (8.5 – 9)	XL (9.5 – 10)
Product Code	02000CHP30-S	02000CHP30-M	02000CHP30-L	02000CHP30-XL

Standards and certifications

ANSI/AAMI PB70 -Level 2 and level 3

EN 13795 - Standard Surgical Gown

EN 388:2016+A1:2018 Clause 6.1, 6.2, 6.3, 6.4

EN ISO 374-1:2016+A1:2018/Type C Level 6

EN ISO 374-2:2019

EN ISO 374-4:2019 Degradation % -9.5

EN ISO 374-5:2016 Protection against Bacterial and Fungi

EN ISO 21420:2020







Back tying belt





Upper tie



Velcro closure




Rubber cuff



100 gowns per box

100,000 gowns per 40 feet container





















State of Israel Ministry of Health משרד הבריאות Amar No. 33640003







Sample Description

Manufacturer





Test Report SL52035260159601TX

Page 1 of 11 Date:November 04, 2020

IMPERTEC INDUSTRIES LTD BARKAN ROAD, BARKAN INDUSTRIAL AREA, 4482100 ISRAEL

The following sample(s) was/were submitted and identified on behalf of the client as:

: (A)Surgical gown

Style No. Composition Sample Color

(A)Non-woven fabric : (A)blue : IMPERTEC INDUSTRIES LTD

Sample Receiving Date Testing Period : Nov 02, 2020 : Nov 02, 2020 - Nov 04, 2020

: Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s). Test Result(s)

Test Performed : Selected test(s) as requested by applicant





काष्ट्राप्त व 国际互认 检测 TESTING CNAS L0599

Page 1 of 4

Test Report SL52035260047201TX Da IMPERTEC INDUSTRIES LTD BARKAN ROAD, BARKAN INDUSTRIAL AREA, 4482100 ISRAEL Date:November 04, 2020

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description : (A)Surgical gown Style No. Composition Sample Color Manufacturer (A)Non-woven fabric (A)Blue IMPERTEC INDUSTRIES LTD

Sample Receiving Date Testing Period

Nov 02, 2020 Nov 02, 2020 - Nov 04, 2020

Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s). Test Result(s)

Test Performed Selected test(s) as requested by applicant

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Gan

Sara Guo (Account Executive)









中国认可 TESTING CNAS L0599

Test Report SL52035260047201TX Da
IMPERTEC INDUSTRIES LTD
BARKAN ROAD, BARKAN INDUSTRIAL AREA, 4482100 ISRAEL Date:November 04 2020 Page 1 of 4

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description : (A)Surgical gown Style No. Composition (A)Non-woven fabric Sample Color (A)Blue IMPERTEC INDUSTRIES LTD Manufacturer

Sample Receiving Date

Nov 02, 2020 Nov 02, 2020 - Nov 04, 2020 Testing Period

Test Result(s) : Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Test Performed : Selected test(s) as requested by applicant

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)







中国队可 国际互认 **CNAS L0599**

Test Report SL52035260047201TX Page 3 of 4 Date: November 04, 2020

Barrier performance of each component and final classification commended

	Impact Penetration Test AATCC 42	Hydrostatic Pressure Test AATCC 127(cmH2O)	Resistance to Bacteriophage Phi-X174 ASTM F 1671	Level	Final classification
Area A (Critical zone-front)	0.4	30.5	/	Level 2	
Area B (Critical zone-sleeve)	0.9	29.9	/	Level 2	
Area C(Front)	0.9	27.3	/	Level 2	Level 2
Area D(Back)	0.6	35.5	/	Level 2	
Seam between areas C&D	0.6	39.5	/	Level 2	
Remark: The barrier performance of all critical zone components, including seams and points of attachments, shall be					

- The barrier performance of all critical zone components, including seams and points of attachments, shall be determined:
 The critical zone comprises at least area A and B. The classification of the surgical gown is based on the lower-performing component of the two;
 The back of the surgical gown (area D) may be nonprotective:
 Seams between protective and nonprotective areas have no barrier requirements;
 Seams between two protective areas are required to have at least the barrier performance of the lower-performing area;
 Level 1: Impact Penetration Test-AATCC 42: s4.5g;
 Level 2: Impact Penetration Test-AATCC 42: s1.0g; Hydrostatic Pressure Test-AATCC 127: ≥20cmH₂O;
 Level 3: Impact Penetration Test-AATCC 42: s1.0g; Hydrostatic Pressure Test-AATCC 127: ≥50cmH₂O;
 Level 4: Resistance to Bacteriophage Phi-X174-ASTM F 1671: Pass.
 Only test the barrier performance of the sample. Label request, construction and other classed in ANSI/AAMI PB70-2012 is not checked.

| Therefore | 1 pp. 1 1 pp.







中国认可 国际互认 检测 TESTING CNAS L0599

Test Report SL52035260047201TX Page 2 of 4

Test Result

Liquid Barrier Performance and Classification of Protective apparel and Drapes Intended for Use in Health Care Facilities
(ANSI/AMIN P870-2012)

Section 4.2.1 Water Resistance: Impact Penetration Test (AATCC 42-2017)

As received					
Weight of blotter gained (g)	1#	2#	3#	4#	5#
Area A (Critical zone-front)	0.3	0.4	0.2	0.1	0.3
Area B (Critical zone-sleeve)	0.9	0.5	0.5	0.8	0.7
Area C(Front)	0.9	0.7	0.9	0.7	0.5
Area D(Back)	0.4	0.3	0.4	0.5	0.6
Seam between areas C&D	0.6	0.4	0.6	0.3	0.3

- Level 1: all critical zone components shall have a blotter weight gain of no more than 4.5grams(g) Level 2: all critical zone components shall have a blotter weight gain of no more than 1.0 grams(g)
- Level 3: all critical zone components shall have a blotter weight gain of no more than 1.0 grams(g)

.2.1 Water Resistance: Hvdrostatic Pressure Test

(AATCC 127-2018; Hydrostatic Head; Rate of increase of water pressure:60mbar/min; temp. of distilled water: 21 °C, fabric face side of water)

As received					
Water Column (cmH ₂ O)	1#	2#	3#	4#	5#
Area A (Critical zone-front)	30.5	38.2	33.8	40.3	41.5
Area B (Critical zone-sleeve)	38.4	29.9	34.3	33.0	35.8
Area C(Front)	33.6	40.5	27.3	35.2	30.9
Area D(Back)	35.5	39.9	37.1	36.1	37.9
Seam between areas C&D	39.5	42.3	42.4	45.6	41.1
Remark:					

Level 2: all critical zone components shall have a hydrostatic resistance of at least 20cmH₂O
 Level 3: all critical zone components shall have a hydrostatic resistance of at least 50cmH₂O



1 (0.21) (0.42) (0.21) (0.42)





CNAS L0599

国际互认

Test Report SL52035260047201TX Page 4 of 4 Date: November 04, 2020



The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.



|Thengs:#the-hu;lockenhapschiq 2000 | (#L25:00000 | (#L25:00000) | |
65-25-8002248669196 | #R. 2000 | | (#L25:0000) | |
66-25-8002248669196 | |







中国认可 国际互认 检测 TESTING CNAS L0599 SL52035260159601TX Date:November 04, 2020 Page 1 of 11

IMPERTEC INDUSTRIES LTD BARKAN ROAD, BARKAN INDUSTRIAL AREA, 4482100 ISRAEL

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description : (A)Surgical gown

: RK-3011A : (A)Non-woven fabric : (A)blue : IMPERTEC INDUSTRIES LTD Style No. Composition Sample Color Manufacturer

: Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

: Selected test(s) as requested by applicant Test Performed

SL52035260159601TX Test Report

Date:November 04, 2020 Page 2 of 11

Comment:

Surgical Clothing and Drapes - Requirements and Test Methods - Part	(A)
1: Surgical Drapes and Gowns (EN 13795-1:2019)	
Clause 4 Bursting strength-Dry	M
Clause 4 Bursting strength-Wet	M
Clause 4 Tensile strength-Dry	M
Clause 4 Tensile strength-Wet	M
Clause 4 Liquid penetration-Water Resistance	M
Clause 4 Cleanliness microbial/Bioburden	M
Clause 4 Microbial Penetration-Dry	M
Clause 4 Microbial penetration-Wet	M
Clause 4 Particle release	M

Remark: M=Meet EN 13795-1:2019 standard performance requirement

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)

Dongjing Liu (Authorized Signatory)





